
Shape-Link Documentation

Release 0.1.1

Paul Müller

Feb 05, 2021

CONTENTS:

1	Getting started	3
1.1	Installation	3
1.2	Citing shapelink	3
2	Command-line interface	5
2.1	shape-link	5
2.1.1	run-plugin	5
2.1.2	run-simulator	6
3	Writing Plug-ins	7
4	Code reference	9
4.1	shapelink.shapelink_plugin	9
4.2	shapelink.util	9
4.3	shapelink.shapein_simulator	10
4.4	shapelink.msg_def	10
5	Changelog	11
5.1	version 0.1.1	11
5.2	version 0.1.0	11
5.3	version 0.0.1	11
6	Indices and tables	13
	Python Module Index	15
	Index	17

This is Shape-Link, a Python library and command-line utility for interfacing with Shape-In, the acquisition software for RT-DC. This is the documentation of Shape-Link version 0.1.1.

GETTING STARTED

1.1 Installation

To install Shape-Link, use one of the following methods:

- **from PyPI:** `pip install shapelink`
- **from sources:** `pip install .`

1.2 Citing shapelink

If you use shapelink in a scientific publication, please cite it with:

Philipp Rosendahl, Paul Müller (2021), Shape-Link version X.X.X: Python library for interfacing with Shape-In [Software]. Available at <https://github.com/ZELLMECHANIK-DRESDEN/shapelink>.

COMMAND-LINE INTERFACE

Shape-Link comes with a command-line-interface (CLI) for running plugins.

2.1 shape-link

```
shape-link [OPTIONS] COMMAND [ARGS]...
```

2.1.1 run-plugin

Run a Shape-Link plugin file

Example usages:

```
# run a plugin
shape-link run-plugin plugins/slp_rolling_mean.py
# run a plugin with a simulator thread (for plugin testing)
shape-link run-plugin -w data.rtdc -f image,deform slp_rolling_mean.py
```

```
shape-link run-plugin [OPTIONS] PATH
```

Options

-w, --with-simulator <with_simulator>

Run the Shape-In simulator in the background using the RT-DC dataset specified (used for testing).

-f, --features <features>

Comma-separated list of features to send by the Shape-In simulator; Defaults to all innate features. A list of valid feature names can be found in the dclab docs (Advanced Usage -> Notation).

Arguments

PATH

Required argument

2.1.2 run-simulator

Run the Shape-In simulator using data from an RT-DC dataset file

Example usage:

```
shape-link run-simulator --features image,deform /path/to/data.rtdc
```

```
shape-link run-simulator [OPTIONS] PATH
```

Options

-f, --features <features>

Comma-separated list of features to send by the Shape-In simulator; Defaults to all innate features. A list of valid feature names can be found in the dclab docs (Advanced Usage -> Notation).

Arguments

PATH

Required argument

WRITING PLUG-INS

A Shape-Link plug-in is a Python script with a class derived from `shapelink.ShapeLinkPlugin` and some additional meta data. Let's have a look at this example plugin which prints the rolling mean of a few scalar features to stdout:

```
1 import shutil
2
3 import numpy as np
4
5 from shapelink import ShapeLinkPlugin
6
7 # We use the terminal width to make sure a line doesn't get cluttered
8 # with prints from a previous line.
9 TERMINAL_WIDTH = shutil.get_terminal_size((80, 20))[0]
10
11
12 class RollingMeansPlugin(ShapeLinkPlugin):
13     """Displays a rolling mean of a few scalar features"""
14     def __init__(self, *args, **kwargs):
15         super(RollingMeansPlugin, self).__init__(*args, **kwargs)
16         self.window_size = 100
17         self.scalar_data = {}
18
19     def after_register(self):
20         print(" Preparing for transmission")
21         for feat in self.registered_data_format.scalars:
22             self.scalar_data[feat] = np.zeros(self.window_size) * np.nan
23
24     def after_transmission(self):
25         print("\n End of transmission\n")
26
27     def handle_event(self, event_data):
28         """Handle a new event"""
29         window_index = event_data.id % self.window_size
30         for ii, feat in enumerate(self.registered_data_format.scalars):
31             self.scalar_data[feat][window_index] = event_data.scalars[ii]
32         # print the first three features to stdout
33         msgs = [" Rolling means: "]
34         num_prints = min(3, len(self.registered_data_format.scalars))
35         for ii in range(num_prints):
36             feat = self.registered_data_format.scalars[ii]
37             msgs.append("{}: {:.3g}".format(feat,
38                                           np.mean(self.scalar_data[feat])))
39         line = " ".join(msgs)
40         if len(line) < TERMINAL_WIDTH:
```

(continues on next page)

(continued from previous page)

```
41         line += " " * (TERMINAL_WIDTH - len(line))
42         print(line, end="\r", flush=True)
43
44         return False
45
46
47 info = {
48     "class": RollingMeansPlugin,
49     "description": "Displays a rolling mean of a few scalar features",
50     "name": "Rolling Means",
51     "version": "0.1.0",
52 }
```

The main action happens in the `handle_event` function - your plugin **must** implement at least this function. The two functions `after_register` and `after_transmission` can be used to set things up (e.g. creation of an additional output file) or to tear things down (e.g. closing that file). Use the `__init__` function for defining additional class properties. The `info` dictionary is required so that the plugin can be run via the [Command-line interface](#).

CODE REFERENCE

4.1 shapelink.shapelink_plugin

Receive data in real-time from a Shape-In instance via zmq

class shapelink.shapelink_plugin.**EventData**

class shapelink.shapelink_plugin.**ShapeLinkPlugin** (*bind_to='tcp://*:6666',* *verbose=False*)

Shape-Link plug-in meta class

Parameters

- **bind_to** (*str*) – IP and port to bind to (where Shape-In runs)
- **verbose** (*bool*) – Set to *True* to see additional debugging information.

after_register ()

Called after registration with Shape-In is complete

after_transmission ()

Called after Shape-In ends data transmission

abstract handle_event (*event_data*: shapelink.shapelink_plugin.EventData) → bool

Abstract method to be overridden by plugins implementations

handle_messages ()

Handle messages from Shape-In

Please don't override this function. Use *ShapeLinkPlugin.handle_event()* for your customized plugins.

4.2 shapelink.util

Utility functions

These functions replicate QDataStreams behavior in C++. In PySide2 QDataStream does not accept array type data.

In C++ an array is serialized by writing: 1) Uint32 number of elements 2) type array elements

more significant bytes are written first. (big-endian)

if numpy “to_bytes” is used the native little-endian format appears

shapelink.util.**qstream_read_array** (*stream*: PySide2.QtCore.QDataStream, *datatype*:
numpy.dtype) → numpy.array

Read array data from a stream with a specified type

```
shapelink.util.qstream_write_array (stream: PySide2.QtCore.QDataStream, array:
                                     numpy.array) → int
    Write array data to a stream with a specified type :param stream: :param array: :return:
```

4.3 shapelink.shapein_simulator

Simulate a Shape-In instance

The communication is based on a simple REQ REP pattern all methods return when the transmission was acknowledged by the peer.

```
class shapelink.shapein_simulator.ShapeInSimulator (destination='tcp://localhost:6666',
                                                    verbose=False)

    register_parameters (scalar_hdf5_names=None, vector_hdf5_names=None, im-
                        age_hdf5_names=None, settings_names=None, settings_values=None)
        Register parameters that are sent to other processes

    send_end_of_transmission ()
        Send end of transmission packet

    send_event (event_id: int, scalar_values: numpy.array, vector_values: List[numpy.array], im-
                age_values: List[numpy.array]) → bool
        Send a single event to the other process

shapelink.shapein_simulator.start_simulator (path, features=None, verbose=1)
    Run a Shape-In simulator using data from an RT-DC dataset
```

4.4 shapelink.msg_def

Definitions for message ids (numeric)

```
shapelink.msg_def.MSG_ID_EOT = -10
    end of transmission

shapelink.msg_def.MSG_ID_EOT_ACK = -11
    end of transmission acknowledged

shapelink.msg_def.MSG_ID_REGISTER = -1
    registration

shapelink.msg_def.MSG_ID_REGISTER_ACK = -2
    registration acknowledged
```

CHANGELOG

List of changes in-between Shape-Link releases.

5.1 version 0.1.1

- docs: add section for writing plugins

5.2 version 0.1.0

- feat: rudimentary command-line interface for running plugins

5.3 version 0.0.1

- initial release

INDICES AND TABLES

- `genindex`
- `modindex`
- `search`

PYTHON MODULE INDEX

S

shapelink.msg_def, [10](#)
shapelink.shapein_simulator, [10](#)
shapelink.shapelink_plugin, [9](#)
shapelink.util, [9](#)

Symbols

--features <features>
 shape-link-run-plugin command line
 option, 5
 shape-link-run-simulator command
 line option, 6
 --with-simulator <with_simulator>
 shape-link-run-plugin command line
 option, 5
 -f
 shape-link-run-plugin command line
 option, 5
 shape-link-run-simulator command
 line option, 6
 -w
 shape-link-run-plugin command line
 option, 5

A

after_register() (*shapelink.shapelink_plugin.ShapeLinkPlugin*
 method), 9
 after_transmission()
 (*shapelink.shapelink_plugin.ShapeLinkPlugin*
 method), 9

E

EventData (*class in shapelink.shapelink_plugin*), 9

H

handle_event() (*shapelink.shapelink_plugin.ShapeLinkPlugin*
 method), 9
 handle_messages()
 (*shapelink.shapelink_plugin.ShapeLinkPlugin*
 method), 9

M

module
 shapelink.msg_def, 10
 shapelink.shapein_simulator, 10
 shapelink.shapelink_plugin, 9
 shapelink.util, 9
 MSG_ID_EOT (*in module shapelink.msg_def*), 10

MSG_ID_EOT_ACK (*in module shapelink.msg_def*), 10
 MSG_ID_REGISTER (*in module shapelink.msg_def*), 10
 MSG_ID_REGISTER_ACK (*in module*
 shapelink.msg_def), 10

P

PATH
 shape-link-run-plugin command line
 option, 6
 shape-link-run-simulator command
 line option, 6

Q

qstream_read_array() (*in module shapelink.util*),
 9
 qstream_write_array() (*in module*
 shapelink.util), 9

R

register_parameters()
 (*shapelink.shapein_simulator.ShapeInSimulator*
 method), 10

S

send_end_of_transmission()
 (*shapelink.shapein_simulator.ShapeInSimulator*
 method), 10
 send_event() (*shapelink.shapein_simulator.ShapeInSimulator*
 method), 10
 shape-link-run-plugin command line
 option
 --features <features>, 5
 --with-simulator <with_simulator>, 5
 -f, 5
 -w, 5
 PATH, 6
 shape-link-run-simulator command line
 option
 --features <features>, 6
 -f, 6
 PATH, 6

ShapeInSimulator (class in
 [*shapelink.shapein_simulator*](#)), 10
shapelink.msg_def
 module, 10
shapelink.shapein_simulator
 module, 10
shapelink.shapelink_plugin
 module, 9
shapelink.util
 module, 9
ShapeLinkPlugin (class in
 [*shapelink.shapelink_plugin*](#)), 9
start_simulator() (in *module*
 [*shapelink.shapein_simulator*](#)), 10